

Ta2 ERASER™

The Professional Strength Solution
For Serious Pigment Problems



Advanced Q-Switched Alexandrite Laser System
for Removal of Tattoos, Pigmented Lesions,
and Other Skin Problems



LIGHT AGE
YOUR AESTHETIC LASER SOURCE

EpiCare™ | Ta2Eraser™ | Q-Clear™

Tel: 732-563-0600
Fax: 732-563-1571
Web: www.lightage.com

Ta2

ERASER™

Ideal for Removing Black, Blue, and Green Tattoos, Refractory Tattoos, Nevi, Age Spots, and More

Specifications

Laser Parameters

Laser Medium:	Alexandrite
Wavelength:	755 nm
Operating Parameters:	Q-switched
Pulse Duration:	75 ± 25 ns
Pulse Energy:	Variable to 800 mJ
Fluence (Energy Density):	Variable to 12 J/cm ²
Spot Sizes:	2 – 6 mm*
Repetition Rate:	Variable to 8 Hz
Cooling:	Closed cycle water-to-air

* 3 mm standard



Physical Parameters

Input Power:	220/240V ac single phase, 15A
Dimensions:	
Height:	43"/109 cm
Width:	19"/48 cm
Length:	28"/71 cm
Delivery System:	Articulated Arm
Control System:	Advanced, Digital



LIGHT AGE
YOUR AESTHETIC LASER SOURCE

EpiCare™ | Ta2Eraser™ | Q-Clear™

